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Valores de leucócitos totais (x10 ⁹ cel/L)	1,0 - 25,0
Neutrófilos (%)	10,0 - 90,0
Linfócitos (%)	5,0 - 85,0
Monócitos (%)	1,0 - 18,0
Eosinófilos (%)	0,0 - 18,0
Basófilos (%)	0,0 - 5,0

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
SÍMBOLOS UTILIZADOS

INSTRUCTIONS FOR USE**INTENDED USE**

The spinit® BC disc is a disposable *in vitro* diagnostic medical device. It is a single-use test, used with the spinit® instrument to perform a quantitative determination of total white blood cells and 5-part differential counts in capillary whole blood samples of human origin.

As an alternative to capillary blood a sample of venous whole-blood can be used. To be used by healthcare professionals.

PRINCIPLES OF THE METHOD/ASSAY

The spinit® BC test is used to determine the total and 5-part differential leucocyte counts from a small blood sample (8µl) collected using the spinit® strip in a 7-minute time period. The microfluidic disc allows for automated sample processing: erythrocyte lysis and white blood cell specific staining. The determination of the leucocyte profile is based on a series of images acquired through an optical microscopy module that captures the specific morphology and the cellular components of each cell type and feeds them into the spinit® system software.

CALIBRATION

The spinit® instrument comes pre-loaded with standard calibration data. All spinit® BC test lots have a unique identification number which are coded together with lot specific calibration data in the product label barcode. A plug-and-play barcode reader is provided to this purpose.

COMPOSITION

Each package contains:

- 20 spinit® BC discs/tests;
- 20 spinit® strips sample collection;
- 1 spinit® BC test IFU (Instruction for Use PL-012-01);
- 1 spinit® strip IFU (Instruction for Use PL-011-01).

Each disc contains the following reagents: saponin, methylene blue, trehalose and an hydrophilic reagent. Necessary materials required but not supplied: sterile lancets approved for capillary sample collection.

WARNING AND PRECAUTIONS

When using capillary blood, the analysis of the disc must start immediately. If the test does not begin within 1 minute, the sample material may dry or coagulate. An information code may be displayed.

- Read carefully the applicable Instructions for Use before performing the test;
- Wear the appropriate PPE (personal protective equipment) while handling samples. Always follow the good laboratory practices;
- Do not use the disc if the pouch has been opened for more than 10 minutes;
- Do not use the test if the pouch or the disc have been damaged;
- Do not use the disc after the expiration date or if the disc have not been stored as recommended;
- Do not re-use the disc;
- Failure to comply with the above recommendations may cause the test to fail and an informative code will be displayed to the user;
- All discs should be disposed of immediately after use. Proper handling and disposal methods should be followed in accordance with the applicable country laws.

STORAGE AND HANDLING

The spinit® BC test disc shelf-life is indicated on the pouch printed label, requiring the discs to be stored in closed pouches at room temperature, consider a dry, clean, well-ventilated area. DO NOT REFRIGERATE OR FREEZE.

SPECIMEN COLLECTION AND PREPARATION

Capillary blood sampling is described in the spinit® strip Instruction for Use PL-011-01.

The capillary blood sample requires no preparation and should be tested immediately.

The EDTA venous whole blood sample should be obtained using analytical methods chosen according to the test requirements, in accordance with good laboratory practices.

The spinit® BC test is ready to use and requires no reagent preparation.

PROCEDURE

Please read spinit® user manual for proper use of the system.
Please read spinit® BC quick guide for detailed information on how to perform a test.
Training videos are available on the instrument's home screen.

1. Remove disc from the box 2. Scan barcode from pouch label and follow on-screen instructions
3. Remove disc from pouch and place it on a flat surface with label facing up 4. Insert the spinit® strip with the sample into the disc. Consult the spinit® strip Instruction for Use PL-011-01 5. Once you have the disc ready, click **Next** on-screen and place the disc on the tray 6. Follow on-screen instructions.

LIMITATIONS OF THE METHOD/ ASSAY

Blood samples with immature leucocytes and nucleated erythrocytes may influence the performance of total and differential WBC count (please refer to section "INTERFERENCE").

Samples with total and differential WBC values outside the test measurement range may impact an increased test failure rate and/or deviating results from reported performance.

PERFORMANCE CHARACTERISTICS**TEST MEASUREMENT RANGE**

The spinit® BC test measurement range for total leucocyte count is 1-25 x 10⁹ cells/L when operated as per recommended procedure. The range for the 5-part differential is presented as a percentage of the total counts.

Total Leucocyte count (x10 ⁹ cells/L)	1.0 - 25.0
Neutrophils (%)	10.0 - 90.0
Lymphocytes (%)	5.0 - 85.0
Monocytes (%)	1.0 - 18.0
Eosinophils (%)	0.0 - 18.0
Basophils (%)	0.0 - 5.0

Table 1 – Quantitative measurement range for spinit® BC test

Results below the Lower Limit of Quantification will appear as lower than (≤) the measurement range defined for the total and 5-part differential leucocyte counts, and similarly results above the Upper Limit of Quantification will appear as higher than (≥) the quantitative range.

METHOD/ASSAY COMPARISON

Results from the method comparison study performed according to CLSI Document EP9-A3 on the spinit® BC instrument and Sysmex XN1000 (standard laboratory reference method) for the total and 5-part differential WBC count.

Parameter	N	SD (x10 ⁹ cel/L)	CV (%)	Bias (%)	R
Leucocytes	196	0.32	6.1	1.3	0.995
Neutrophils	196	0.53	7.1	2.8	0.984
Lymphocytes	196	0.33	9.9	-0.4	0.971
Monocytes	196	0.12	14.7	0.7	-
Eosinophils	196	0.05	-	1.2	-

Table 2 – Method comparison analytical performance

PRECISION**Inter-day**

Precision studies were conducted according to the CLSI EP5-A3 guideline. The results were obtained from duplicate measurements of control solutions distributed over the analytical range, using 2 spinit® instruments and 1 production lot of spinit® BC discs over a period of 14 days (n = 60). See table 3.

Sample	Parameter	Density (x10 ⁹ cel/L)	SD (x10 ⁹ cel/L)	CV (%)	Repeatability (%)
1	Leucocytes	19.7	1.27	6.5	5.9
	Neutrophils	10.2	0.74	7.2	6.2
	Lymphocytes	6.5	0.46	7.0	6.9
	Monocytes	2.4	0.18	7.3	5.9
	Eosinophils	0.5	0.07	13.9	12.7
2	Leucocytes	10.4	0.63	6.0	5.6
	Neutrophils	5.7	0.40	7.0	6.1
	Lymphocytes	3.6	0.20	5.5	4.8
	Monocytes	0.8	0.09	12.4	10.8
	Eosinophils	0.3	0.03	12.7	12.4

Table 3 – Precision of spinit® BC test

STANDARDIZATION/ACCURACY

A range of clinical samples were tested replicated 20 times on the same day with the same lot of spinit® BC test and 1 spinit® instrument. A summary of the results is presented in table 4.

Parameter	Density (x10 ⁹ cel/L)	SD (x10 ⁹ cel/L)	CV (%)	Bias (%)
Leucocytes	6.2	0.40	6.5	0.6
Neutrophils	3.2	0.23	7.2	2.1
Lymphocytes	2.1	0.22	10.8	-0.6
Monocytes	0.5	0.08	16.8	-1.6
Eosinophils	0.3	0.07	22.2	6.7

Table 4 – Accuracy of spinit® BC test

INTERFERENCE

The spinit® system informs the user for the presence of immature leucocytes and nucleated red blood cells with specific sample flags. Refer to spinit® instrument user manual for information on sample flags. For specific performance characteristics, please refer to www.biosurfit.com or contact your supplier for more information.

CONTROL MATERIAL

Each spinit® BC test disc contains multiple integrated quality control verification steps to ensure reliability of results. The control kit for spinit® BC is recommended for routine quality control testing. Consult the Instruction for Use PL-019-01 of spinit® BC Control kit (ref. 910201). Use QC mode for control testing. See spinit® User Manual PL-006-05 for further details.

ISO13485 certified company.

SYMBOLS USED: