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APLICAÇÃO

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INSTRUCTIONS FOR USE

INTENDED USE

spinit® strips are manufactured to be used together with spinit® products. spinit® strips are designed for the collection of different types of samples (capillary blood, venous blood, serum, plasma, controls, calibrators) and transfer to spinit® discs. The spinit® strips are to be used by properly trained healthcare professionals only in accordance with these instructions.

spinit® strip is a sampling method specifically designed and developed for spinit® tests, to be used in a near patient test spinit® system.

PRODUCT DESCRIPTION


spinit® strips are non-sterile low sample volume collection strip with a predefined nominal volume for achieving correct amount of sample needed to a specific in vitro spinit® test.

The spinit® strip is comprised of an elongated body having a handling portion disposed at one end and arranged to provide for gripping and handling means and a sample collecting portion disposed at the other end and arranged to provide for means to receive a volume of a liquid sample.

The spinit® strip is used to collect a volume of sample (7.5 - 8 µL) by touching the blood drop or liquid.

The device further comprises a frangible portion disposed between the handling portion and the sample collecting portion, wherein the frangible portion is configured to enable complete separation of the sample collecting portion from the sample transfer device.

spinit® STRIP SAMPLE COLLECTION

spinit® STRIP HANDLING

Please read these instructions, your spinit® test instructions and your spinit® instrument User Manual before using spinit® strips.

STEPS

- Open the vial and remove a collection strip; close the vial immediately after (Vial contains a bag with desiccant);
- Prepare the sampling area by cleaning the finger using disinfectant.
- Firmly prick the patient's finger and wipe away the first drop of blood;

- Collect the second drop of blood by bringing the tip of the strip to the surface of the blood drop as shown in figure.
- Observe the sample area is fully loaded. If not, it is possible to refill;
- Clean the strip by gently wiping each side with tissue paper;

- After cleaning the strip, observe if sample area is filled, if not, discard the strip and collect a fresh sample using a new collection strip;
- Insert the strip by sliding the strip tip into the disc sample well;
- Fold back to break the strip and leave the sample inside the disc.

IMPORTANT!

- Do not clean the strip in a vertical position as this will remove the sample;
- Do not use tissue paper on the tip of the strip;
- Capillary blood samples must be collected from one finger prick (do not squeeze finger for multiple collections from same finger prick);
- Capillary blood samples must be analysed immediately after collection.

STORAGE AND HANDLING

Store the strips in their original vial in a cool, dry place, away from sunlight, at room temperature 2 to 30° C (35,6° - 86° F). Do not freeze.

Note: If delivered within a spinit® test box, the spinit® strip can be stored together with the respective box at the specific temperature indicated in each spinit® test IFU/label. When you take a strip from the vial, close the cap immediately. Use the strip immediately.

Do not transfer the spinit® strip to a new vial or any other container.

Check the expiry date on the vial. If the expiry date has passed do not use the strips.

Note: All expiration dates are printed in Year/Month/Day format.

MATERIALS PROVIDED WITH 1 (ONE) VIAL

- 20 spinit® strip
- 1 silica bag

LIMITATIONS

Use each spinit® strip once. Do not reuse. Do not bend, cut or damage the spinit® strip. spinit® strip should not be used with other samples than the ones indicated.

PERFORMANCE CHARACTERISTICS

The performance of spinit® strip have been evaluated with the spinit® disc and spinit® instrument.

COMPOSITION

Made of PET material.

SAFE DISPOSAL OF THE DEVICE

Disposal in accordance with local legislation.

CUSTOMER SERVICE

Adverse events will not occur if the indications on the current Instruction for use are followed, and the product is handled by a professional.

Additional information, as well as instructions for use of spinit® products, can be found or request at:

All suspected serious adverse events must be reported to the distributor, manufacturer or regulatory authority in the Member State where the user is located.

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